

# Octo-Pos Thorax Positioner



- + COMFORT
- + MODULARITY
- + FLEXIBILITY

# Octo-Pos

## Modular and Extensible



The inside has a 45° bevel



Octo-Pos is made of foam, coated with water-repellent single-component, high strength, breathable, formaldehyde-free, biologically degradable polyurethane resin, it does not contain any toxic additives


**Max comfort**  
thanks to a new one  
rounded design

The new design of the Octo-Pos positioner offers rounded shapes, especially on the critical contact points. This way, pressure damage is avoided, particularly in the long-term interventions.

Ideal also for pediatric patients

## Total flexibility thanks to the modular additional spacers

The positioner is equipped with standard removable 5cm spacers that allow the optimal setting in order to position the patient. For further sizing, 2.5 cm spacer kits are available.



## Octo-Posis:

### Versatile


Only one device for all patients

### Extensible

The flexibility of the device allows it to be extensible, obtaining different sizes.

### Modular

Each module is covered with Velcro and this allows maximum positioning flexibility.



**Ideal** also for non  
conventional postures

Ideal for interventions that require to expose not  
reachable anatomical areas in conventional postures.


# Octo-Pos

## Maximum stability and modularity

Octo-Pos allows maximum stability thanks to the base covered with Velcro®. This solution allows personalized modules positioning and guarantees stable anchoring avoiding detachments and slipping.



Large Velcro® base



Octo-Pos is equipped with a set of spacers, in order to guarantee a perfect adaptation of the device to the patient.

### BioMatrix

Via Maffini, 23 - 29010 San Giuliano di Castelvetro (PC)  
Tel. 0523.826619 Fax. 0523826728  
e-mail: [info@biomatrix.it](mailto:info@biomatrix.it)

All rights reserved. Changes may be made without prior notice.  
BioMatrix S.r.l. is not to be considered responsible for damages  
resulting from the lack or inaccuracy of the information reported here.

[www.biomatrix.it](http://www.biomatrix.it)

